

# TMS Background and Motivation


**1 in 8**

Canadians suffer from depression

**20% – 30%** of cases are **treatment resistant**



**TMS** (Transcranial Magnetic Stimulation) is a noninvasive treatment for treatment-resistant depression




## How does it work?

- The TMS coil generates a **magnetic field** at a specific hotspot (location of stimulation on the cerebral cortex)
- The magnetic field induces an **emf** which drives **electric signals**
- Signals then modulate **neural activity** and influence **brain function**

However, the current method where the coil is manually moved by a doctor results in a **lack of precision & efficiency**

## Our **Solution**

Based on Dr. Takemi's design, we aim to build an autonomous cobot system that can move and maintain its position with **higher precision and efficiency** than a human operator.




[1] [www.researchgate.net/figure/The-schematic-of-ECoG-electrode-in-the-presence-of-TMS-coil\\_fig1\\_233330682](https://www.researchgate.net/figure/The-schematic-of-ECoG-electrode-in-the-presence-of-TMS-coil_fig1_233330682)

[2] <https://soterixmedical.com/research/hd/neuronavigation>

[3] <https://www.universal-robots.com/products/ur3-robot/>

# System Overview

- - - - Dynamic connections      ——— Static connections



## IR Reflecting Head Accessory

Allows IR camera to track patient's head in real time

## IR Camera

Tracks coil and patient location

## IR Reflecting Orbs


Allows IR camera to track TMS coil in real time

## UR3e Cobot

Cobot used in this project

## TMS Coil

Delivers TMS to patient





# Calibration and Registration

## 1) Cobot Calibration & Patient Registration



## 2) TMS Treatment







First, we calibrate the cobot by registering points in the patient's **MRI anatomical space** to the **physical space** seen through the IR camera.


With this, we are able to autonomously position the coil to any point on the MRI using the cobot and real time IR data.

During this process, users will be able to monitor the cobot's progress through a Graphical User Interface (**GUI**).

This will display real time **IR coordinates** of the Cobot, Coil, and Patient, overlaid upon the existing MRI.



Nasion




We can calibrate the cobot by mapping points between anatomical space and image space

[1] <https://www.fieldtriptoolbox.org/faq/coordsys/>

[2] <https://pharmrev.aspetjournals.org/content/pharmrev/74/4/918.full.pdf>

# Implementation



## Next Steps:

- Iterate on Dr. Takemi's system with a combined force and IR camera cobot controller
- Manufacture 3D printed TMS coil holder, robot fixing stand, and subject pointer
- Refine and improve cobot force control loop
- Integrate all system components to create a useable platform for psychiatric researchers in the NINET Lab

[1] <https://optitrack.com/cameras/v120-duo/>

[2] <https://optitrack.com/software/motive/>

[3] <https://medium.com/analytics-vidhya/getting-started-with-ros-overview-installation-and-ros-computational-graph-model-e94d7a16187f>

[4] <https://www.universal-robots.com/products/ur3-robot/>